

a measurable difference...

Sole Adhesion Tester

Model: S0010

This well proven instrument was designed to measure the strength of the adhesion of stuck-on soles at the toe and heel in the shoe factory, but is equally useful in the laboratory. A steel plate is used as a cantilever to pull the sole of the footwear away from the upper. The maximum force reading of the dial indicator is 100kg.

Application:

- Footwear

Features:

- Analogue Gauge
Range 0 - 100kg x 1

Standards:

- AS/NZ 2210.2
- BS 5131
- ISO 782/1

Options:

- Heel attachment to suit
120mm wide heels

Connections:

- N/A

Dimensions:

- **H:** 100mm
- **W:** 90mm
- **D:** 300mm
- **Weight:** 6kg



10-11 Colrado Court, Hallam, Victoria 3803 Australia.

Phone: +61 3 9708 6885 Fax: +61 3 9708 6770

Website: www.idminstruments.com.au Email: ids@idminstruments.com.au

All sizes and dimensions shown on this data sheet are averages only and images are as accurate as the medium allows.
Products and details shown may change without notice as IDM Instruments reserves the right to make alterations to its product range from time to time.